



## TECHNICAL BULLETIN : PRODUCT INFORMATION

The following contains a detailed list of the compatibility of specific NuAire basestands with the Ergotron accessory.

### NU-500 Series BSC's

The NuAire 500 Series of Biological Safety Cabinets are Ergotron compatible with the following basestands.

The following adjustable (hydraulic) basestands are compatible with the NU-500-142 Ergotron Workstation and NU-M432 Ergotron Laptop Tray.

#### All are compatible

-300	31" with Leg Levelers 115V	NU-500-125
-300	31" with Leg Levelers 230V	NU-500-133
-300	31" with Casters 115V	NU-500-129
-300	31" with Casters 230V	NU-500-137
-400	31" with Leg Levelers 115V	NU-500-126
-400	31" with Leg Levelers 230V	NU-500-134
-400	31" with Casters 115V	NU-500-130
-400	31" with Casters 230V	NU-500-138
-500	31" with Leg Levelers 115V	NU-500-127
-500	31" with Leg Levelers 230V	NU-500-135
-500	31" with Casters 115V	NU-500-131
-500	31" with Casters 230V	NU-500-139
-600	31" with Leg Levelers 115V	NU-500-128
-600	31" with Leg Levelers 230V	NU-500-136
-600	31" with Casters 115V	NU-500-132
-600	31" with Casters 230V	NU-500-140

The following telescoping basestands are compatible with the NU-500-141 Ergotron Workstation and NU-M429 Ergotron Laptop Tray.

**Must be shipped after 1/2017**

-300	31" with Leg Levelers	NU-500-100
-300	31" with Casters	NU-500-108
-300	31" BSSC with Leg Levelers	NU-500-104
-300	31" BSSC with Casters	NU-500-112
-400	31" with Leg Levelers	NU-500-101
-400	31" with Casters	NU-500-109
-400	31" BSSC with Leg Levelers	NU-500-105
-400	31" BSSC with Casters	NU-500-113
-500	31" with Leg Levelers	NU-500-102
-500	31" with Casters	NU-500-110
-500	31" BSSC with Leg Levelers	NU-500-106
-500	31" BSSC with Casters	NU-500-114
-600	31" with Leg Levelers	NU-500-103
-600	31" with Casters	NU-500-111
-600	31" BSSC with Leg Levelers	NU-500-107
-600	31" BSSC with Casters	NU-500-115

## NU-400 Series BSC's

The NuAire 400 Series of Biological Safety Cabinets are Ergotron compatible with the following basestands.

The following telescoping basestands are compatible with NU-500-141.

### Must be shipped after 3/2017

-300	28" with Leg Levelers	NU-400-211
-300	28" with Casters	NU-400-219
-300	28" BSSC with Leg Levelers	NU-400-215
-300	28" BSSC with Casters	NU-400-223
-400	28" with Leg Levelers	NU-400-212
-400	28" with Casters	NU-400-220
-400	28" BSSC with Leg Levelers	NU-400-216
-400	28" BSSC with Casters	NU-400-224
-500	28" with Leg Levelers	NU-400-213
-500	28" with Casters	NU-400-221
-500	28" BSSC with Leg Levelers	NU-400-217
-500	28" BSSC with Casters	NU-400-225
-600	28" with Leg Levelers	NU-400-214
-600	28" with Casters	NU-400-222
-600	28" BSSC with Leg Levelers	NU-400-218
-600	28" BSSC with Casters	NU-400-226

## NU-140/201/240 Laminar Airflow Cabinets

The NuAire 140, 201, and 240 Laminar Airflow Cabinets are Ergotron compatible with the following basestands.

The following telescoping basestands are compatible with NU-200-141.

### Must be shipped after 3/2017

-300	28" with Leg Levelers	NU-200-133
-300	28" with Casters	NU-200-137
-300	28" OPM with Leg Levelers	NU-200-143
-400	28" with Leg Levelers	NU-200-134
-400	28" with Casters	NU-200-138
-400	28" OPM with Leg Levelers	NU-200-144
-500	28" with Leg Levelers	NU-200-135
-500	28" with Casters	NU-200-139
-500	28" OPM with Leg Levelers	NU-200-145
-600	28" with Leg Levelers	NU-200-136
-600	28" with Casters	NU-200-140
-600	28" OPM with Leg Levelers	NU-200-146

The following motorized basestands are compatible with NU-200-163.

**Must be shipped after 3/2017**

-300	28" with Leg Levelers 115V	NU-200-147
-300	28" with Casters 115V	NU-200-149
-300	28" with Leg Levelers 230V	NU-200-148
-300	28" with Casters 230V	NU-200-150
-400	28" with Leg Levelers 115V	NU-200-151
-400	28" with Casters 115V	NU-200-153
-400	28" with Leg Levelers 230V	NU-200-152
-400	28" with Casters 230V	NU-200-154
-500	28" with Leg Levelers 115V	NU-200-155
-500	28" with Casters 115V	NU-200-157
-500	28" with Leg Levelers 230V	NU-200-156
-500	28" with Casters 230V	NU-200-158
-600	28" with Leg Levelers 115V	NU-200-159
-600	28" with Casters 115V	NU-200-161
-600	28" with Leg Levelers 230V	NU-200-160
-600	28" with Casters 230V	NU-200-162

## **NU-340 Laminar Airflow Cabinets**

The NU-340-100 Ergotron Workstation is compatible with all NU-340 Laminar Flow Cabinets. However, unlike other Ergotron kits the NU-340 version is not just a bolt-on accessory and requires additional holes to be drilled. Please see drawings CD-000073, CD-000088, & CD-000156 for further details.

## NU-813 Containment Cabinets

The NuAire 813 Containment Cabinets are Ergotron compatible with the following basestands.

The following telescoping basestands are compatible with NU-800-146.

### Must be shipped after 3/2017

-300	28" with Leg Levelers	NU-800-128
-300	28" with Casters	NU-800-131
-300	28" Extended Height with Leg Levelers	NU-800-140
-300	28" Extended Height with Casters	NU-800-143
-400	28" with Leg Levelers	NU-800-129
-400	28" with Casters	NU-800-132
-400	28" Extended Height with Leg Levelers	NU-800-141
-400	28" Extended Height with Casters	NU-800-144
-600	28" with Leg Levelers	NU-800-130
-600	28" with Casters	NU-800-133
-600	28" Extended Height with Leg Levelers	NU-800-142
-600	28" Extended Height with Casters	NU-800-145

The following motorized basestands are compatible with NU-800-147.

### Must be shipped after 3/2017

-300	MOTORIZED BS W/CSTR 115V 3'	NU-800-148
-300	MOTORIZED BS W/CSTR 230V 3'	NU-800-149
-300	MOTORIZED BS W/LL 115V 3'	NU-800-150
-300	MOTORIZED BS W/LL 230V 3'	NU-800-151
-400	MOTORIZED BS W/CSTR 115V 4'	NU-800-152
-400	MOTORIZED BS W/CSTR 230V 4'	NU-800-153
-400	MOTORIZED BS W/LL 115V 4'	NU-800-154
-400	MOTORIZED BS W/LL 230V 4'	NU-800-155
-600	MOTORIZED BS W/CSTR 115V 6'	NU-800-156
-600	MOTORIZED BS W/CSTR 230V 6'	NU-800-157
-600	MOTORIZED BS W/LL 115V 6'	NU-800-158
-600	MOTORIZED BS W/LL 230V 6'	NU-800-159