


BASE STAND ASSEMBLY



A

$\frac{3}{8}$ "-16 x 2 $\frac{3}{4}$ " Bolt
 Leg Bushing
 $\frac{3}{8}$ " Washer

$\frac{3}{8}$ " Washer
 $\frac{3}{8}$ " Lock Washer
 $\frac{3}{8}$ "-16 Acorn Nut

B

$\frac{1}{4}$ "-20 x $\frac{5}{8}$ " Bolt
 $\frac{1}{4}$ " Lock Washer
 $\frac{1}{4}$ " Washer

$\frac{1}{4}$ " Washer
 $\frac{1}{4}$ " Lock Washer
 $\frac{1}{4}$ "-20 Acorn Nut

C

$\frac{1}{4}$ "-20 x $\frac{5}{8}$ " Bolt
 $\frac{1}{4}$ " Lock Washer
 $\frac{1}{4}$ " Washer



NuAire, Inc. | nuaire.com